# Nasojejunal & nasogastric feeding tubes

### Description:

These feeding tubes are designed for continuous or intermittent nasogastric or nasointestinal enteral feeding, for neonatal, paediatric and adult patients. The tubes are used to feed, flush and deliver medicines and to aspirate enteral contents to confirm correct position. The tubes are equipped with a male ENFit connector compliant with UNI CEI EN ISO 80369-3 European Standard, and are made of biocompatible clear polyurethane, they are fully radio-opaque and visible under X-Ray. Tubes are graduated every centimeter to ease insertion and control possible displacement; these marks are only approximate. NGP and LG tubes have a closed atraumatic tip, while NJP tubes have an open tip which allows them to be passed over an endoscopically or radiologically placed guidewire if required. Selected lengths and sizes are available in longest weighted versions; all weighted tubes have a closed tip. NGP and NJP ranges are equipped with a stainless steel guidewire.

Tubes are suitable for long-term use (over 30 days) and have been tested for use up to 90 days. Care should be taken with ventilated patients and local hospital policies should be followed.

Using an appropriate ENFit syringe, the tube should be flushed before and after giving feeds and administering medications. The external portion of the tube should be cleaned and dried daily. The devices must be used exclusively by healthcare professionals.

**ENTERAL USE ONLY**  
**CAREFULLY READ THE INSTRUCTIONS BEFORE USE.**

**Manufacturer**  
**HMC Premedical S.p.A.**  
Via Tonino Morandi 16,  
41037 Mirandola (MO) Italy  
tel.: +39 0535 22704  
fax.: +39 0535 609546  
**Distributor**  
**Medicina Ltd.**  
Units 1-4,  
Rivington View Business Park,  
Station Road, Bolton  
BL6 5BN UK

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