# ENFit® nasogastric feeding tube

## Instructions For Use

**DEVICE DESCRIPTION:**  
The ENFit Nasogastric Feeding Tube equipped with Male ENFit Connector in line with EN ISO 80369-3 European Standard.

These tubes are designed for nasogastric or nasointestinal, continuous or intermittent enteral feeding. Tubes are made of polyurethane or silicone, they are biocompatible and radiopaque, or have a radiopaque line to be easily located by X-rays. The tubes are graduated to ease insertion and control any possible displacement.

The long-term tubes NGP range have guide wires and the short-term tubes CG & RT ranges do not.

**INDICATIONS AND INTENDED USE:**  
The short-term tubes (CG) are mainly used in neonates and children but can be used in adults who need short-term enteral support.

Short-term wide bore tubes (RT) are used to aspirate and feed adults and children. Typically they are used in intensive care, high dependency units, theatres, crash trolleys and surgical wards.

Short-term tubes should not be used for longer than 30 days. Long-term tubes (NG) are polyurethane tubes with guidewires. They are used for long term support to feed adults and children. Long term tubes should not be used for longer than 60 days

When necessary the device should be replaced.

**CONTRAINDICATIONS:**  
Do not use Nasogastric Tubes in the case of:  
(1) Gastrointestinal failure,  
(2) Intestinal obstruction,  
(3) Acute abdomen,  
(4) Active gastrointestinal bleeding.

Do not use CG or NGP/C range in the case of:  
1) High volume decompression in emergency situations  
2) Decompression for high volume or high viscous fluid

Clinical judgement and/or a risk assessment based on the patient's needs should always be performed by a healthcare professional when choosing to use an RTx/C tube under the following circumstances:  
1) High volume decompression in emergency situations  
2) Decompression for high volume or high viscous fluid

**INTENDED USER:**  
The device must be used exclusively by healthcare professionals, or suitably trained personnel.

**STORAGE CONDITIONS**  
Keep Dry  
Keep away from Sunlight

**DISPOSAL:**  
For the disposal of the Nasogastric Tubes follow the local Hospital/Healthcare facility protocols or as directed by the Healthcare professional.

### ENTERAL USE ONLY CAREFULLY READ THE INSTRUCTIONS BEFORE USE.

**CAUTION: FEDERAL LAW (USA)** restricts this device to sale by or on the order of a physician. Follow Physicians instructions for use of this device.

|   |  |   |   |   |
|---|--|---|---|---|
|  |  |  |  |  |
| Single Use  | Sterilized using Ethylene Oxide  | Do not resterilize  | Caution: This device is for sale by or on order of a physician                      | Consult Instructions for use  |
|  |  |  |  |  |
| Do not use with phthalates, phthalates or BPA                                     | Does not contain Natural Rubber Latex  | Keep dry  | Keep away from sunlight   | Do not use if packaging is damaged  |
|  |  |   |   |   |

**Medicina Ltd.**  
Units 2-4, Rivington View Business Park,  
Station Road, Bolton BL6 5BN UK  
Tel: 01204 695050

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**Cleaning protocol:**

Clean the work area, wear protective gloves, use a new ENFit Connector 50 or 60ml enteral syringe, flush with warm water before and after aspiration/delivery, dry the outside of the feeding tube, check the integrity of the tube for signs of damage/deterioration.

For further cleaning instructions for non-health care professionals, please refer to the GEDSA website.
**http://stayconnected.org/wp-content/uploads/2018/05/ENFit-Cleaning-Procedures-All-Tubes.pdf**

**DIRECTIONS:**

**Insertion**

The tube is stretched to remove any memory retained from the packaging, which prevents the guide-wire sticking. The length of the tube is estimated by measuring from the sternum to the tip of nose added to the tip of nose to ear lobe. This length should be noted against the scale on the tube. For insertion, the tip of the tube is lubricated with water and gently inserted into the nose with the head bent slightly forward. If possible, the passage of the tube can be aided by the patient swallowing a sip of water. If the tube needs to be re inserted into the patient within its intended duration of use the cleaning protocol should be repeated.

**Confirmation of position**

Correct position is confirmed by aspiration of stomach contents using an appropriate pH strip or by X-ray. The clear short-term tubes have a radiopaque stripe and the long-term tubes are fully radiopaque. The indicator is used by attaching a ENFit enteral syringe (20 or 50mls) into the guide wire holder ENFit Connector or the purple ENFit connector of the tube. The connections should be secure and airtight. A small quantity of fluid should be aspirated and its pH tested. When pH testing you should refer to pH strip manufacturer’s IFU.

**Securing the tube**

Having carefully removed the guide wire the position of the tube should be noted. Then the tube can be secured to the nose using an appropriate Nasogastric Plaster. In addition the tube can be secured to the face under the cheekbone with tape.

**Connections**

Feeding and medicines connections are made to the purple ENFit connector on the tube. Aspiration or suction equipment can be connected to the funnel of the RT range of tubes. All tubes should be flushed before and after medicines, after feeds and every four hours, with 20mls of water in adults and 10-20mls of water in children and 2-5mls in neonates. (Care should be taken with fluid restricted patients). If the tube blocks a small quantity of fizzy water can be used but the guidewire should never be re-introduced whilst the tube is still in the patient. When securing the cap to the tube connector, only use a ¼ (quarter) turn to effect the seal of the tube.

An Enteral Feeding Pump to control the flow of enteral nutrition can also be used with feeding tubes.

|   |
|---|
| <b>WARNINGS/ PRECAUTIONS:</b>   |
| <b>Caution</b> NOTHING should be introduced down the tube before gastric placement has been confirmed. DO NOT FLUSH the tube before gastric placement has been confirmed. The internal guide wire should NOT be lubricated before gastric placement has been confirmed. |
| <b>Caution</b> Feeding should not start if there is any doubt about the position of the tube. The indicator strip should show a definite colour change. When pH testing you should refer to the pH strip manufacturer’s IFU and your own organisation’s procedures.     |
| <b>Caution</b> Small syringes can generate high pressures sufficient to burst the tube. Syringes may also cause the tube to collapse if excessive pressures are generated. Aspiration should not be continued if resistance is felt.                                    |
| <b>Caution</b> If the strap attaching the cap to the connector breaks, firmly screw the cap into the connector and replace the tube as soon as possible.  |
| <b>Caution</b> Do not use tools in conjunction with connection or disconnection of the connector or the cap.  |
| <b>Caution CG or NGP/C range only.</b> Do not use this device for high volume decompression in emergency situations, those at high risk of aspiration and/ or for decompression for high volume or high viscous fluid.  |
| <b>Caution RT range.</b> Do not use the ENFit connector on this device for aspiration, except where drawing up aspirate for confirmation of tube placement.   |

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