

## **TECHNICAL DATA SHEET**

Name	ENFit Home Sy	ringe						
Product Code	LHExx (Non-St	erile)						
Shelf Life	5 Years							
Lot Number forma	t HMYYxxx (Y is	HMYYxxx (Y is for year, xxx is a sequential batch number)						
Expiry and								
Manufactured date	e YYYY/MM	YYYY/MM						
format								
	General	All components (product, Bag, and boxes) should be free of damage, no tears or holes.						
	Blister	Product is provided in a single blister pack.						
	Вох	White box with print shall be free of any major damage.						
	Outer Carton	Shipping carton shall be free of any major damage.						
	Syringe	Foreign matter under 30001x – 7001x Illuminance. Free from foreign particles and matters.						
		Syringe free from defects; fins, burrs, plastic flow and material shortage.						
Appearance		Syringe barrel shall be transparent to see fiducial line clearly. Syringe Barrel shall be						
		legible with clear and even lines.						
		All print should be legible and clearly seen.						
		Plunger will be attached to black piston and shall be intact.						
		Shall be free from	n blooming and foreigr	n matters. The surface shall	be free from cracks,			
	Piston	cutting, damage and uneven surfaces.						
		Outer diameter of the piston shall be equal or shorter than the outer diameter of the						
		piston sealing ring.						
Labelling and	Label Information	shall match the copies in artwork master file.						
Symbols	Symbols	Shall match the copies in artwork master file.						
Material	Barrel, Plunger & Tip							
Composition	Piston	Silicone Rubber (Latex-Free)						
	Lubricant	Silicone Oil.						
Packaging		Multiple package (see table below). QML (minimum delivery quantity): 1 Box						
		PRODUCT CODE	CAPACITY	BOX CONTENT	CARTON CONTENT			
	Bag	LHE01LD	1ml (low-dose)	100	2400			
		LHE25LD	2.5ml (low-dose)	100	2000			
		LHE05	5ml	100	1600			
		LHE10	10ml	100	1000			
		LHE20	20ml	80	800			
		LHE30	30ml	80	640			
		LHE60	60ml	60	240			
		LHE100	100ml	30	120			



		pH Value difference	e between extra	ict and blan	k solution shall not exceed	1.0 pH		
			ution: Purified v			·		
	pH Value		tion: Purified wa		•			
Chemical		Note, the quantity	Note, the quantity of extract solution is based on the nominal capacity of the syringe.					
Performance	Readily	The syringe extrac	The syringe extract solution and equal volume of blank solution to indirect titrations as					
	Oxidisable	per Standard GB/	KM <sub>n</sub> O <sub>4</sub> consumption sha	ll not exceed				
	substance	0.5ml.						
	IR Spectra	The spectra of the syringe materials correspond to the reference spectra.						
	Conical	No evident shaking or movement between standard gauge and conical fitting as per ISO						
	Fitting	594.						
		There shall be no air leakage when syringe is under negative pressure 88Kpa for 60s + 5s						
	Leakage	at the contacting position between barrel and piston. The piston will not detach from						
		plunger.						
		Nominal capacity	Tolerance on graduated ca			-Maximum Dead Space ml		
		of syringe, V ml	Less than half nominal		Equal to or greater than			
			capacity	( . <b>f</b>	half nominal capacity			
Physical Performance		V < 2	±(1.5% of V + 2% of expelled volume)		± 5% of expelled volume	0.07		
		2 ≤ V < 5	±(1.5% of V + 2% of expelled volume)		± 5% of expelled Volume	0.07		
	Capacity Tolerance	5 ≤ V < 10	±(1.5% of V +1% of expelled volume)		± 4% of expelled Volume	0.075		
	lorenance	10 ≤ V < 20	±(1.5% of V +1% of expelled volume)		± 4% of expelled Volume	0.10		
		20 ≤ V < 30	±(1.5% of V +1% of expelled volume)		± 4% of expelled Volume	0.15		
		30 ≤ V < 50	±(1.5% of V +1% of expelled volume)		± 4% of expelled Volume	0.17		
		50 ≤ V	±(1.5% of V +1% of expelled volume)		± 4% of expelled Volume	0.20		
		Nominal Capacity of syringe, V ml	Initial Force, Fs Max. N	Mean Force, F Max N	Maximum Force F Max N	Minimum Force, F min N		
Sliding Performance		V< 2	10	5	(2,0 measured ) or (measured+1,5N), whichever is higher	(0,5 measured ) or (measured– 1,5N), whichever is the lower		
		2 ≤ V < 50	25	10	(2,0 measured ) or (measured +1,5N), whichever is higher	(0,5 measured ) or (measured– 1,5N), whichever is the lower		
		50 ≤ V	30	15	(2,0 measured ) or (measured +1,5N), whichever is higher	(0,5 measured ) or (measured– 1,5N), whichever is the lower		



	When syringe is filled with water up to its nominal capacity and kept vertically, the				
	plunger shall not be moved by its own weight.				
Product Drawing					
	Latex	This product is free from Latex			
	Phthalates	This product is free from Phthalates			
	Animal Derived Substances	This product is free from Animal Derived Substances			
Product Inclusions	Bis Phenol A	This product is free from Bis Phenol A			
	Asbestos	This product is free from Asbestos			
	PCB's	This product is free from PCB's			
Storage	To be stored at a minimum of of 80%.	To be stored at a minimum of 5°C and a maximum of 35°C with a maximum relative humidity of 80%.			
Disposal	The user must follow the leg	The user must follow the legal regulations regarding disposal of hospital waste.			
	Medicina Itd Units 1-4 Rivington View Business park				
Manufacturer	Station Road, Blackrod Bolton BL6 5BN UK				
Classification	Class IIa Annex V & VII	Class IIa Annex V & VII			
Conformity Assessment Rou	te: MDD 93/42/EEC				