

## **TECHNICAL DATA SHEET**

Name	ENFit Enteral S	Syringe						
Product Code	LPExx (Sterile)							
Shelf Life	5 Years							
Lot Number forma	at HMYYxxx (Y is	HMYYxxx (Y is for year, xxx is a sequential batch number)						
Expiry and Manufactured dat format								
Appearance	General	All components (product, Bag, and boxes) should be free of damage, no tears or holes.						
	Blister	Product is provided in a single blister pack.						
	Вох	White box with print shall be free of any major damage.						
	Outer Carton	Shipping carton shall be free of any major damage.						
		Foreign matter under 30001x – 7001x Illuminance. Free from foreign particles and matters. Syringe free from defects; fins, burrs, plastic flow and material shortage.						
	Syringe	Syringe barrel shall be transparent to see fiducial line clearly. Syringe Barrel shall be legible with clear and even lines. All print should be legible and clearly seen.						
		Plunger will be attached to black piston and shall be intact.						
	Piston	<ul><li>Shall be free from blooming and foreign matters. The surface shall be free from cracks, cutting, damage and uneven surfaces.</li><li>Outer diameter of the piston shall be equal or shorter than the outer diameter of the</li></ul>						
		piston sealing ring.						
Labelling and Symbols	Label Information	Shall match the copies in artwork master file.						
	Symbols	Shall match the copies in artwork master file.						
Material	Barrel, Plunger & Tip	& Tip Polypropylene.						
Composition	Piston	Silicone Rubber (Latex-Free)						
	Lubricant	Silicone Oil.						
		Multiple package (see table below). QML (minimum delivery quantity): 1 Box						
		PRODUCT CODE	CAPACITY	BOX CONTENT	CARTON CONTENT			
		LPE01LD	1ml (low-dose)	100	2400			
Packaging		LPE25LD	2.5ml (low-dose)	100	2000			
	Bag	LPE05	5ml	100	1600			
		LPE10	10ml	100	1000			
		LPE20	20ml	80	800			
		LPE30	30ml	80	640			
		LPE60	60ml	60	240			
		LPE100	100ml	30	120			
	Charillin							
Biological	Sterility	Shall be sterile.						



		pH Value difference	pH Value difference between extract and blank solution shall not exceed 1.0 pH						
	PH Value		• Extract solution: Purified water + Sample.						
		Blank Solution: Purified water.							
Chemical									
Performance									
	Oxidisable	per Standard GB/T14233.1, the amount of KM <sub>n</sub> O₄ consumption shall not exceed							
	substance	0.5ml.							
	IR Spectra		The spectra of the syringe materials correspond to the reference spectra.						
Physical Performance	Conical Fitting	No evident shaking or movement between standard gauge and conical fitting as per ISO 594.							
	Fitting	There shall be no air leakage when syringe is under negative pressure 88Kpa for 60s + 5s							
	Leakage	at the contacting position between barrel and piston. The piston will not detach from plunger.							
		Nominal capacity	Tolerance on gr	aduated ca	pacity				
		of syringe, V ml	Less than half r capacity		Equal to or greater than half nominal capacity	-Maximum Dead Space ml			
		V < 2	±(1.5% of V + 2% of expelled volume)		± 5% of expelled volume	0.07			
		2 ≤ V < 5	±(1.5% of V + 2% of expelled volume)			0.07			
	Capacity Tolerance	5 ≤ V < 10	volume)		± 4% of expelled Volume	0.075			
		10 ≤ V < 20	±(1.5% of V +1% of expelled volume)		± 4% of expelled Volume	0.10			
		20 ≤ V < 30	±(1.5% of V +1% of expelled volume)		± 4% of expelled Volume	0.15			
		30 ≤ V < 50	±(1.5% of V +1% of expelled volume)		± 4% of expelled Volume	0.17			
		50 ≤ V	±(1.5% of V +1% of expelled volume)		± 4% of expelled Volume	0.20			
Sliding Performance		Nominal Capacity of syringe, V ml	Initial Force, Fs Max. N	Mean Force, F Max N	Maximum Force F Max N	Minimum Force, F min N			
		V< 2	10	5	(2,0 measured ) or (measured+1,5N), whichever is higher	(0,5 measured ) or (measured– 1,5N), whichever is the lower			
		2 ≤ V < 50	25	10	(2,0 measured ) or (measured +1,5N), whichever is higher	(0,5 measured ) or (measured– 1,5N), whichever is the lower			
		50 ≤ V	30	15	(2,0 measured ) or (measured +1,5N), whichever is higher	(0,5 measured ) or (measured– 1,5N), whichever is the lower			



	When syringe is filled with water up to its nominal capacity and kept vertically, the			
	plunger shall not be moved by its own weight.			
Product Drawing				
	Latex	This product is free from Latex		
	Phthalates	This product is free from Phthalates		
Product Inclusions	Animal Derived Substances	This product is free from Animal Derived Substances		
	Bis Phenol A	This product is free from Bis Phenol A		
	Asbestos	This product is free from Asbestos		
	PCB's	This product is free from PCB's		
Storage	To be stored at a minimum of 5°C and a maximum of 35°C with a maximum relative humidity of 80%.			
Disposal	The user must follow the legal regulations regarding disposal of hospital waste.			
Manufacturer	Medicina Itd Units 1-4 Rivington View Business park Station Road, Blackrod Bolton BL6 5BN UK			
Classification	Class IIa Annex V & VII			
Conformity Assessment Route:	MDD 93/42/EEC			