

SUGRBAG®

COOL MED WALLET Activation



THANK FOR PURCHASING SUGRBAG® COOL MED WALLET !

These instructions are only a guide. Duo to varied climatic and environmental conditions, SUGRBAG® COOL MED wallets will respond differently according to location. If activated and used sensibly, they will continue to respond to your needs over a long period of time. With age, the length of immersion time will increase. However, as long as the inner pouch's size expand, the wallet will function properly. Otherwise , it will need to be replaced.

Activation Instructions

1. Immerse the inner pouch, in cold water.
2. Soak the wallet for 5-8 minutes
3. Spread the crystals evenly through the pouch. The activated pouch should remain pliable.
4. Towel the inner pouch and place medication before puting it inside the outer cover.

Please Notes:

- The SUGRBAG[®] COOL MED inner pouch should be re-immersed in water when the size of the inner becomes smaller and hot (in continuous use, succeeding immersion periods will be shorter, i.e. 3-4 minutes). It is extremely important that the correct amount of water is absorbed during immersion because over-soaking causes the wallet to become too swollen to be used! In his case, the wallet must be allowed to dry naturally until it returns to a usable size.
- The SUGRBAG[®] s need to breathe. They work by evaporation and need to be in contact with the air. The outer cover allows this to happen enhancing the stability of the pouch's temperature. Carrying in hand luggage or in a pocket is satisfactory.
-
- Care should be taken when towel drying your SUGRBAG[®] COOL MED pouch after immersion in water, as there might be some slight color loss from the fabric of brightly-colored pouches.
- The SUGRBAG[®] COOL MED with medication, can be placed in a refrigerator but **DO NOT PLACES IN A FREEZER.**

DO NOT CARRY OR STORE IN AN AIRTIGHT OR WATERPROOF CONTAINER WHEN ACTIVATED!

Storage Instructions

When not in use, remove the inner pouch from the outer cover and allow it to completely dry out naturally until the gel returns to crystal form. To keep the crystals separate, give the pouch an occasional shake whilst drying out. When completely dry, your SUGRBAG[®] COOL MED should be stored in dry conditions until next required. The drying out period can take 2-4 weeks depending on the environmental climate.

To accelerate the drying out period, various methods can be used:

1. place the inner pouch near a heat source such as a heat vent
2. hang the inner pouch on a clothes line if the weather is suitable
3. leave the inner pouch on a window sill in a sunny position

DO NOT PLACE THE SUGRBAG[®] COOL MED IN A TUMBLE DRYER!

Cleaning Instructions

The outer cover of the wallet can be washed in warm water using a very mild detergent. The inner wallet can be spot cleaned by using a soft bristle brush. Dip the brush into slightly soapy water and spot clean the wallet. You can then rinse the wallet under running water without damaging the crystals inside. Alternatively, you can use a non-bleach anti-bacterial cleaning spray or an anti-bacterial kitchen or baby wipe to clean your wallet.

DO NOT PLACE THE SUGRBAG[®] COOL MED WALLET IN A WASHING MACHINE!