


# TECHNICAL DATA SHEET

## ENFit Infant short term Feeding Tube

**Reference:** CGxx/xxx

**Trade Name:** ENFit Infant short term Feeding Tube Polyurethane



### **PRODUCT DESCRIPTION:**

ENFit Infant short term Feeding Tube Polyurethane, available from 4FR to 14FR

### **INTENDED USE:**

The device must be used exclusively by healthcare professionals, or suitably trained personnel. The Short Term Tubes (CG) are mainly used in neonates and paediatrics but can be used with adults who need short term enteral support.

### **RAW MATERIAL:**

Tube: Polyurethane and barium sulphate

Connector: ABS

Strap: Polyethylene

Cap : Polypropylene

### **PACKAGING:**

Primary packaging: Single product in Tyvek pouch.

Secondary Packaging (Box): White Cardboard box – quantities refer to table 1

Tertiary Packaging (Shipping Carton): Brown corrugated cardboard carton– quantities refer to table 1

**Table 1 - Product Range:**

	Outer Diameter	Inner Diameter	Length of Tube	Box quantity	Carton Quantity
CG4/40	1.38	0.6	400(+10/-2)	50	500
CG4/50	1.38	0.6	500(+10/-2)	50	500
CG4/50W	1.38	0.6	500(+10/-2)	50	500
CG4/80	1.38	0.6	810(+10/-2)	50	500
CG4/90	1.38	0.6	910(+10/-2)	50	500
CG5/40	1.7	0.9	400(+10/-2)	50	500
CG5/50	1.7	0.9	500(+10/-2)	50	500
CG5/80	1.7	0.9	810(+10/-2)	50	500
CG6/40	2	1.3	400(+10/-2)	50	500
CG6/50	2	1.3	500(+10/-2)	50	500
CG6/50W	2	1.3	500(+10/-2)	50	500
CG6/80	2	1.3	810(+10/-2)	50	500
CG6/80W	2	1.3	810(+10/-2)	50	500
CG6/100	2	1.3	1100(+10/-2)	50	500
CG6/120	2	1.3	1210(+10/-2)	50	500
CG8/40	2.7	1.5	400(+10/-2)	50	500
CG8/50	2.7	1.5	500(+10/-2)	50	500
CG8/80	2.7	1.5	810(+10/-2)	50	500
CG8/80W	2.7	1.5	810(+10/-2)	50	500
CG8/120	2.7	1.5	1210(+10/-2)	50	500
CG10/80	3.3	2.15	810(+10/-2)	25	500
CG10/80W	3.3	2.15	810(+10/-2)	25	500
CG10/120	3.3	2.15	1210(+10/-2)	25	500
CG12/80	4	2.8	810(+10/-2)	25	500
CG12/80W	4	2.8	810(+10/-2)	25	500
CG12/120	4	2.8	1210(+10/-2)	25	500
CG14/80	4.7	3.2	810(+10/-2)	25	500

**STERILISATION:**

EtO Sterilisation is undertaken in accordance with the requirements of the current regulations  
The device should not be re-sterilised.

**SHELF LIFE:**

Shelf life 5 years from sterilisation

**USAGE:**

Single use device

**STORAGE:**

Standard storage procedures and conditions.

**DISPOSAL:**

The user must follow the legal regulations regarding disposal of hospital waste.

**MANUFACTURER:**

Medicina Ltd  
Units 2-4 Rivington View Business Park, Station Road, Blackrod, BL6 5BN

**EC REPRESENTATIVE:**

HMC Premedical S.p.A.  
Via Bosco 1/3, 41037, Mirandola (MO) – Italy

**NOTIFIED BODIES:**
**SGS Belgium NV – CE1639**

Certificate Number: GB19/964577

**SGS United Kingdom – UKCA0120**

Certificate number – GB22/00000566

**CLASSIFICATION:**

Class IIa

**CONFORMITY ASSESSMENT ROUTE:**

MDD 93/42/EEC Annex V

UK MDR 2002 (as amended) Part II

**REVISION HISTORY:**

Revision No:	Date:	Summary of changes	Issued by:
01	07/05/2020	First issue with CE1639	Joanne Schmidt
02	17/02/2021	Addition of EC Representative	Joanne Schmidt
04	21/12/2022	Include UKCA and update packaging	Joanne Schmidt