


HMC PREMEDICAL S.p.A.
VIA TONINO MORANDI, 16
41037 MIRANDOLA (MO) ITALY
TEL. 0535 22704 - FAX. 0535 609546
INFO@HMCGROUP.IT

(€ 0051

TECHNICAL DATASHEET

ENFit Enteral blunt fill needle with filter



CODE	BOX pcs.	CARTON pcs.	
MN02L	50		



HMC PREMEDICAL S.p.A.
VIA TONINO MORANDI, 16
41037 MIRANDOLA (MO) ITALY
TEL. 0535 22704 - FAX. 0535 609546
INFO@HMCGROUP.IT

CE

PRODUCT COMPOSITION:

- ENFit small bore male connector in purple ABS, produced in compliance with UNI CEI EN ISO 80369-3 norm.
- 2. Blunt needle 18G x 38 mm in Stainless Steel and Polypropylene, endowed with a 5 µm fluid filter (not shown).

The device has a needle with **blunt tip at 45 degrees**; compared to a common hypodermic needle, it requires a force only 2.5 times higher to perforate a rubber membrane, and 10 times higher to pierce the human skin. The needle is also endowed with a **5 µm fluid filter**, which enhances patient safety by reducing glass particles contamination of fluids drawn from glass ampoules.

INTENDED USE:

Aspiration of fluids from glass ampules.

PACKAGING:

Single package in medical grade paper/ Film PP-PE Blister. Multiple package: boxes of 50 pcs; cartons of 200 pcs.

PRODUCTION PROCESS:

The device is manufactured in accordance with HMC Premedical S.p.A. Quality System and in compliance with the requirements of the standard EN ISO 13485.

CONTROL ON THE PRODUCT:

In all stages of processing, according to internal procedures and to sampling plans defined by norm ISO 2859-1.

CLASSIFICATION:

Class IIa sterile.

STERILIZATION:

Product sterilized by ethylene oxide according to a validated course of treatment in accordance with current regulations.

Shelf life 5 years from the date of sterilization.

For single-use only.

Non re-sterilizable.

STORAGE:

Standard storage procedures and conditions.

DISPOSAL:

The user must follow the legal regulations regarding disposal of hospital waste.

WARNING:

The device must be used exclusively by healthcare professionals.

REGISTRATION TO ITALIAN REPERTOIRE:

CND: **A0880** RDM: **1572532**

REV. Rev.	MODIFICHE Changes	DATA Date	EMESSO Issued	VERIFICATO e APPROVATO Verified and Approved
00	First issue – English language	01/08/2017	E. Benassi	D. Bosetti
01	Data sheet update	13/10/2017	E. Benassi	D. Bosetti
02	Product composition update	25/08/2018	E. Benassi	D. Bosetti
03	Packaging update; layout revision	20/01/2020	E. Benassi	S. Tralli
				1