

Qufora IrriSedo Mini

For low volume transanal irrigation

Description

Qufora IrriSedo Mini includes:

- **A.** Single-use cone with hydrophilic coating that is activated when dipped in water
- B. Pump

Accessories

Single-use cones

Extension tube

Features

- The soft single-use cone is coated so when dipped in water the surface becomes lubricated and it is ready to use straight away
- The cone is single-packed in a pouch that can be used to attach and detach the cone hygienically, which can also be used as a disposal bag for the used cone
- · The procedure can be performed using only one hand
- The manual pump has a non-return valve that prevents used water from running back in the pump
- If a little extra water is needed, or using both hands is prefered, a single-use extention tube is available

Shelf life

5 years from production

Disposal

Qufora IrriSedo Mini does not require special disposal protocols and may be disposed of in household waste according to local guidelines and recommendations

Classification and sterilization

Qufora IrriSedo Mini is a Class I Device, non-sterile

Regulations and standards:

- MDR (Regulation 2017/745) CE
- UK MDR 2002 UKCA
- ISO 13485:2016 (certified)

Qufora IrriSedo Mini



Hydrophilic coated cone

Single-use cone is ready to use when dipped in water



Pump

Contains 90 ml water



Extension tube

Accessory item to deliver an increased volume of water (approx 120ml), and can aid dexterity

Material description

Parts	Material	Softener	
Cone			
Clear cone/coated	TPE	No softener	
Coating	PVP	No softener	
Pump			
Pump	PVC, phthalate-free	DINCH	
White NRV valve out and in	Hard PVC	No softener	
Mesh for valve	Hard PVC	No softener	
Black valve balls	Nitrile	No softener	
Tube	PVC, phthalate-free	DEHT/DOTP	
Purple connector	Hard PVC	No softener	
Extension tube			
Extension tube	PVC, phthalate-free	DEHT/DOTP	
Purple connector	Hard PVC	No softener	

Dimensions	Length	Diameter	Volume
Cone	75 mm	Tip is Ø8 mm	
Pump	170 mm	59 mm	90 ml
Extension tube	60 cm	9,2 mm	

For more information, please visit www.qufora.com