


## TECHNICAL SHEET

<b>Reference</b>	SF01
<b>Trade Name:</b>	Surfactant Kit
<b>Product Description</b>	Surfactant Kit with needle
<b>Intended Use</b>	For administration of pulmonary Surfactant
<b>Raw Materials:</b>	<ol style="list-style-type: none"> <li>1. M.L.L Superlock Connector - blue Rigid – PVC. DEHP Free.</li> <li>2. Graduated Tube Ø 1x1.7 mm (CH5) Length 215 mm in Polyurethane with blue radiopaque stripe.</li> <li>3. Syringe with Superlock F.L.L Connector in PP.</li> <li>4. M.L.L Superlock connector in blue Rigid – PVC. DEHP Free.</li> <li>5. M.L.L Connector – ABS.</li> <li>6. Steel needle 18G x 38mm with rounded-out endpoint and protective capsule.</li> </ol>
	
<b>Packaging</b>	Single package in paper PP-PE film blister. Multiples packed in boxes of 25
<b>Shelf Life:</b>	5 years from the date of sterilisation
<b>Usage</b>	Single use device
<b>Sterilisation</b>	EtO Sterilisation of the Surfactant kit is undertaken in accordance with the requirements of the current regulations
	Device should not be re-sterilised
<b>Storage</b>	Normal storage procedures and conditions

<b>Disposal</b>	The user must follow the legal regulations regarding disposal of hospital waste.
<b>Manufacturer</b>	HMC Premedical Spa Via Tonino Morandi 16 41037 Mirandola Italy
<b>Classification</b>	Class IIa Sterile
<b>Conformity Assessment Route:</b>	MDD 93/42/EEC, Annex VII
<b>Notified Body</b>	IMQ NB Number: 0051
<b>Certificate Ref</b>	1668/MDD