

TECHNICAL SHEET

Reference	SF02
Trade Name:	Needle free Surfactant Kit
Product Description	Surfactant Kit without needle
Characteristics	
Intended Use	For administration of pulmonary surfactant
Raw Materials:	 M.L.L Superlock connector - blue Rigid – PVC. DEHP free. Graduated Tube Ø 1x1.7 mm (CH5) Length 215 mm in Polyurethane with blue radiopaque stripe. Syringe with superlock F.L.L Connector in PP. M.L.L Superlock connector in blue Rigid – PVC. DEHP free. Strap – PVC. F.L.L Superlock cap – PVC. DEHP free. M.L.L. connector – ABS. Vial cap Ø 13mm -PC
	Sippedox 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1
Packaging	Single package in paper PP-PE film blister.
Shalf Life.	Multiples packed in boxes of 25
Shelf Life:	5 years from the date of sterilisation
Usage	Single use Device
Sterilisation	EtO Sterilisation of the Surfactant kit is undertaken in accordance with the requirements of the current regulations
_	Device should not be Re-sterilised
Storage	Normal storage procedures and conditions



Disposal	The user must follow the legal regulations regarding disposal
	of hospital waste.
<mark>IFU</mark>	
Manufacturer	HMC Premedical Spa
	Via Tonino Morandi 16
	41037 Mirandola
	Italy
Classification	Class IIa Sterile
Conformity Assessment Route:	MDD 93/42/EEC, Annex VII
Notified Body	IMQ NB Number: 0051
Certificate Ref	1668/MDD