

WELLAND MEDICAL

Leaders in Stomacare Technology

FINISHED PRODUCT DATA SHEET

Product: **Welland Adhesive Remover**

Supplied as:	<u>Format</u>	<u>Code</u>
	Wipe	WAD 050
	Aerosol	WAS 040
	Pump Dispenser	WAS 050

Description:

Welland Adhesive remover is a non sting solution for skin application that is designed to easily remove adhesive from the skin leaving no residue. This will then enable re-sticking on the same surface a few minutes after removal.

The key ingredients are a silicone solvent, hexamethyldisiloxane, a "wetting out agent" and, in the aerosol version, a propellant.

The solution is designed to trap between the adhesive and the skin, releasing the adhesive as it moves.

Technical Specification:

Active constituents	-	Siloxanes
Flash point	-	-6°C (DN3.3)
Ignition Point	-	352°C
Refractive Index	-	1.3790
Colour	-	Colourless
Odour	-	Characteristic Odour

November 2006

Page 1 of 2

\\wmnt1\data\Private\QA\QA Data\SPCS\FINISHED PRODUCT DATA SHEETS\Adhesive Remover.doc

Welland Medical Limited
7/8 Brunel Centre, Newton Road, Crawley, West Sussex RH10 9TU

Tel: +44 (0) 1293 615455/ 615422 Fax: +44 (0) 1293 615411/ 615412 Email: mail@wellandmedical.com Web: www.wellandmedical.com



INVESTOR IN PEOPLE

A CliniMed Group Company

Registered Office: Cavell House, Knaves Beech Way, Loudwater, High Wycombe, Bucks. HP10 9QY Registered in England 2170535 VAT No. 475 6939 87

WELLAND MEDICAL

Leaders in Stomacare Technology

Wipe Format:

Dose rate	-	1.9ml (1.5g) per sachet
Wipe size	-	70mm x 150mm
Sachet material	-	Paper/Polyethylene/Aluminium/Surlyn
Sachet Size	-	60mm x 50mm (x 2 perforated)
Boxed	-	50 per box

Aerosol Format:

Contents of can	-	50 ml
Welland Adhesive Remover	-	40.3ml
Propellant	-	9.7ml
Can Material	-	Shell – Aluminium Cap – Polypropylene Moving parts – Nylon

Pump Dispenser Format:

Contents of Can	-	50ml
Welland Adhesive Remover	-	50ml
Can material	-	Shell – Aluminium Cap – Polypropylene

November 2006

\\wmnt1\data\Private\QA\QA Data\SPCS\FINISHED PRODUCT DATA SHEETS\Adhesive Remover.doc

Welland Medical Limited
7/8 Brunel Centre, Newton Road, Crawley, West Sussex RH10 9TU

Tel: +44 (0) 1293 615455/ 615422 Fax: +44 (0) 1293 615411/ 615412 Email: mail@wellandmedical.com Web: www.wellandmedical.com



INVESTOR IN PEOPLE

A CliniMed Group Company

Registered Office: Cavell House, Knaves Beech Way, Loudwater, High Wycombe, Bucks. HP10 9QY Registered in England 2170535 VAT No. 475 6939 87