


TECHNICAL SHEET

Reference	IN01/02
Trade Name:	pH Indicator strips
Product Description	pH Indicator strips for human gastric aspirate testing
Intended Use	The pH Indicator strips are used to determine the correct placement/misplacement of nasogastric feeding tubes by testing the pH value of human aspirate.
Raw Materials:	1. Pure paper (Cellulose) 2. Indicator dye



Packaging	Multiple package (see table below). QML (minimum delivery quantity): 1 Box Reference Features: Unit <table border="1" style="width: 100%;"> <thead> <tr> <th>code</th> <th>Box/Tube</th> </tr> </thead> <tbody> <tr> <td>IN01</td> <td>200 per box</td> </tr> <tr> <td>IN02</td> <td>100 per tube</td> </tr> </tbody> </table>	code	Box/Tube	IN01	200 per box	IN02	100 per tube
code	Box/Tube						
IN01	200 per box						
IN02	100 per tube						
Shelf Life:	Non sterile						
Usage	Single use device						
Storage	Normal storage procedures and conditions						
Disposal	The user must follow the legal regulations regarding disposal of hospital waste.						
Manufacturer	Medicina ltd Units 1-4 Rivington View Business park Station Road Blackrod						

	Bolton BL6 5BN
Classification	Class I
Conformity Assessment Route:	MDD 93/42/EEC, Annex VII