


TECHNICAL SHEET

Reference	Oral tip syringe																														
Trade Name:	Oral tip syringe																														
Product Description	Oral tip syringes are used to deliver medication and feed orally. Used by patients in the hospital and community environment.																														
Intended Use	To deliver feed/Medication orally																														
Raw Materials:	<ul style="list-style-type: none"> • COMPONENT 1: Barrel-Polypropylene • COMPONENT 2: Plunger-Polypropylene • COMPONENT 3: Sealing ring-Silicon rubber • COMPONENT 4: Lubricant-Medical grade silicone oil • COMPONENT 5: Oral tip - Polypropylene 																														
																															
Packaging	Multiple package (see table below). QML (minimum delivery quantity): 1 Box <table border="1" style="margin-top: 10px; width: 100%;"> <thead> <tr> <th style="background-color: #d9ead3;">CODICE</th> <th style="background-color: #d9ead3;">ML</th> <th style="background-color: #d9ead3;">PACKAGE Pcs.</th> <th style="background-color: #d9ead3;">BOX Pcs.</th> </tr> </thead> <tbody> <tr> <td>OT005</td> <td>0.5ml</td> <td>1</td> <td>100</td> </tr> <tr> <td>OT01</td> <td>1ml</td> <td>1</td> <td>100</td> </tr> <tr> <td>OT05</td> <td>5ml</td> <td>1</td> <td>100</td> </tr> <tr> <td>OT10</td> <td>10ml</td> <td>1</td> <td>100</td> </tr> <tr> <td>OT20</td> <td>20ml</td> <td>1</td> <td>80</td> </tr> <tr> <td>OT25</td> <td>2.5ml</td> <td>1</td> <td>100</td> </tr> </tbody> </table>			CODICE	ML	PACKAGE Pcs.	BOX Pcs.	OT005	0.5ml	1	100	OT01	1ml	1	100	OT05	5ml	1	100	OT10	10ml	1	100	OT20	20ml	1	80	OT25	2.5ml	1	100
CODICE	ML	PACKAGE Pcs.	BOX Pcs.																												
OT005	0.5ml	1	100																												
OT01	1ml	1	100																												
OT05	5ml	1	100																												
OT10	10ml	1	100																												
OT20	20ml	1	80																												
OT25	2.5ml	1	100																												
Shelf Life:	5 years from the date of sterilization																														
Usage	Single use device																														
Storage	Normal storage procedures and conditions																														
Disposal	The user must follow the legal regulations regarding disposal of hospital waste.																														
Manufacturer	Medicina ltd Units 1-4 Rivington View Business park Station Road Blackrod Bolton BL6 5BN																														
Classification	Class IIa Annex V & VII																														
Conformity Assessment Route:	MDD 93/42/EEC																														