





*Feel the difference.
Make a difference.*



mediven[®] cotton

mediven cotton is the first medical compression stocking with organic cotton from medi. The combination with functional fibres provides a thermoregulating effect. The natural fibre cotton is considered particularly soft – and therefore suitable for sensitive and allergic skin which has been dermatologically tested and given top marks. The use of organic cotton and recycled polyamide yarn provides additional benefits for the environment and the people involved.

mediven cotton. Naturally well cared for.

mediven[®] cotton

Naturally well cared for.

medi tip¹:

Due to the material thickness, the fit and the product characteristics, we recommend mediven cotton for:

- Women / men
- Normal weight / overweight
- Light to heavy tendency to venous oedema
- Also for soft connective tissue

Indications²

- Chronic venous disease stages I to II according to Widmer or COs – C4 according to CEAP
- Initial phase following varicose therapy
- Congestion problems and oedema during pregnancy
- Postoperative reperfusion oedema
- Post-traumatic, postoperative, occupational, medicinally induced oedema
- Oedema due to immobility
- Superficial venous thrombosis
- Condition after thrombosis, post-thrombotic syndrome
- Thrombosis prophylaxis in mobile patients
- Inflammatory dermatoses of the legs
- Nausea and dizziness in pregnancy

Features

- First medical compression stocking with organic cotton – lined on the inside of the stocking for direct skin contact
- Polyamide yarn with recycled content for effective and resource-saving compression




- Thermoregulating effect – cooling in hot temperatures and warming in cold weather
- Ideal compression solution for patients with atopic and sensitive skin – dermatologically proven via dermatest
- Opaque look

¹ Please note that this recommendation does not replace a doctor's diagnosis or professional advice from a specialist medical retailer.

² Selecting the ideal compression stockings depends on the indication and the individual needs of the patients.

Washing instructions

Easy-care 


Suitable for washing machine and drier.



Double woven with organic cotton yarn



Organic cotton – directly on the skin



Open toe – very narrow border

Product benefits



Material


Medium

Compression class

CCL 1, 2

Standard colours

S = Standard



* Delivery period up to 10 working days.

mediven cotton

open toe

Styles	AD	AD pte.	AG / NoB	AG pte./NoB	AT	AT pte.
	calf	calf petite	thigh with topband	thigh petite with topband	panty	panty petite
Class 1	Y540x01-7	Y547x01-7	Y569x01-7	Y559x01-7	Y510x01-7	Y517x01-7
Class 2	Y640x01-7	Y647x01-7	Y669x01-7	Y659x01-7	Y610x01-7	Y617x01-7

closed toe

Styles	AD	AD pte.	AG / NoB	AG pte./NoB	AT	AT pte.
	calf	calf petite	thigh with topband	thigh petite with topband	panty	panty petite
Class 1	Y540x11-7	Y547x11-7	Y569x11-7	Y559x11-7	Y510x11-7	Y517x11-7
Class 2	Y640x11-7	Y647x11-7	Y669x11-7	Y659x11-7	Y610x11-7	Y617x11-7


See the last page for the article numbers

Standard styles		
Standard	Size I – VII, Length standard and petite, topband standard and wide	
Made-to-measure	not available	
Special models	not available	
Border	AD	tricot border
	AG	Silicone dot topband, Silicone micro dot topband sensitive
All topbands are standard styles.		
Gusset	AT	Standard gusset
Foot part	open toe, closed toe, Shoe size 36 – 43	

Material composition			
	Elastane (EL)	Polyamide (PA)*	Organic cotton
CCL 1	51%	22%	27 %
CCL 2	46%	27%	27 %

*With a recycled content of 60 percent

Measurement chart

Standard fit without seam		I	II	III	IV	V	VI	VII
AT gusset	cT	from 60–125	from 60–125	from 60–125	from 60–125	from 60–125	from 60–125	from 60–125
	cH	from 80–135	from 80–135	from 80–135	from 80–135	from 80–135	from 80–135	from 80–135
AG with topband, normal	cG	43–48	45–52	49–56	53–60	56–64	60–68	64–72
AG with topband, wide	cG	49–57	53–62	57–67	61–72	65–77	69–81	73–85
 <p>ℓ=Length c=Circumference</p>	cG	43–57	45–62	49–67	53–72	56–77	60–81	64–85
	cF	39–52	41–56	44–60	47–65	50–69	53–73	57–78
	cE	30–37	33–40	35–43	37–45	39–48	41–51	42–53
	cD	27–33	29–36	32–39	34–42	36–45	38–48	40–50
	cC	28–34	30–37	33–40	35–43	37–46	39–49	41–51
	cB1	23–27	24–29	26–32	29–35	31–37	33–39	35–41
	cB	18–20	20–22	22–24	24–26	26–28	28–30	30–32
	cY	26–31	28–33	29–35	31–37	32–38	33–40	34–42
	cA	17–22	19–24	21–26	23–29	25–32	27–34	28–35
	Length measurements	ℓG regular	72–83 cm	Foot		Highly elastic universal foot size, shoe size 36–43		
	ℓG petite	62–71 cm						
	ℓD regular	39–44 cm						
	ℓD petite	34–38 cm						

medi tip:



medi Butler donning made easier¹ **medi skin and textile care offer support**

More information on the medi accessories from page 124

¹ In general, the following applies for donning and doffing aids: Rabe E et al. Medical compression therapy of the extremities with medical compression stockings (MCS), phlebological compression bandages (PCB), and medical adaptive compression systems (MAC): S2k guideline of the German Phlebology Society (DGP) in cooperation with the following professional associations: DDG, DGA, DGG, GDL, DGL, BVP. Hautarzt 2021;72(2):137-152.