


TECHNICAL SHEET

Reference	NPxx
Trade Name:	Nasogastric Plasters
Product Description	Adult / Paediatric Nasogastric plasters
Intended Use	Nasogastric Plasters are used to secure a nasogastric tube to the nose of a patient. It is as an adhesive plaster with wings that are attached across the bridge of the nose this is a single-use device.
Raw Materials:	<ol style="list-style-type: none"> 1. Spunlaced non-woven fabric 2. Solvent based pressure-sensitive acrylic adhesive 3. Glassine release paper
 <p style="margin-left: 200px;"> NP01 – Adult plaster NP02 – Paediatric plaster </p>	
Packaging	NP01/NP02 Plasters Bag: 100 single plasters Shipping carton 10 bags, 1000 units QML (minimum delivery quantity): 100 units (1 bag) Description of reference: Reference Features: Unit
Shelf Life:	3 years from the date of manufacture
Usage	Single use device
Sterilisation	Non Sterile
Storage	Normal storage procedures and conditions
Disposal	The user must follow the legal regulations regarding disposal of hospital waste.
Manufacturer	Medicina Ltd. Units 1-4, Rivington View Business Park, Station Road, Blackrod, Bolton, BL6 5BN
Classification	Class I Non Sterile