




Juzo Compression Wrap Fodsegment

Art. 6000 LF		Seriestørrelser · Omfang i cm				
			2 / S	3 / M	4 / L	5 / XL
	cB		15,0 – 23,0	20,0 – 28,0	25,0 – 33,0	30,0 – 40,0
	cA'		19,0 – 24,0	24,0 – 29,0	29,0 – 34,0	34,0 – 40,0
	PA	normal	16	18	19,5	19,5
		lang	19	21	22,5	23
		* Cirka 6 cm proximalt fra lilletåens grundled				


Juzo Compression Wrap Underbensegment

Art. 6000 LC/LCM		Seriestørrelser · Omfang i cm						
Længde i cm			2 / S	3 / M	4 / L	5 / XL	6 / XXL	7 / XXXL
	PD							
	AD							
	middel < 40,0 lang > 40,0							
	cC		29,0 – 35,0	34,0 – 40,0	37,0 – 44,0	41,0 – 48,0	46,0 – 53,0	50,0 – 60,0
	max		35,0 – 41,0	40,0 – 46,0	44,0 – 51,0	48,0 – 55,0	53,0 – 60,0	60,0 – 70,0
	cB		18,0 – 21,0	21,0 – 24,0	24,0 – 27,0	27,0 – 31,0	31,0 – 35,0	36,0 – 47,0

Juzo Compression Wrap Knæsegment

Art. 6000 LK		Seriestørrelser · Omfang i cm					
		2 / S	3 / M	4 / L	5 / XL	6 / XXL	7 / XXXL
	cE	30,0 – 36,0	36,0 – 42,0	42,0 – 48,0	48,0 – 55,0	55,0 – 62,0	62,0 – 70,0

Juzo Compression Wrap Lårsegment

Art. 6000 LT		Seriestørrelser · Omfang i cm				
Længde i cm		2 / S	3 / M	4 / L	5 / XL	6 / XXL
	PFG					
	cG	43,0 – 70,0	54,0 – 81,0	65,0 – 92,0	75,0 – 102,0	85,0 – 112,0
	cF	35,0 – 55,0	48,0 – 70,0	60,0 – 82,0	70,0 – 92,0	80,0 – 102,0
		middel < 25,0				
		lang > 30,0				